

1. Part Number: A336-1000

2. Photograph:



3. Product Information:

Description:	Open PharmaScoop® 1000ml
Material of Construction:	316L Stainless Steel
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	1000ml
Overall Length:	395mm
Body Length:	270mm
Body Diameter:	135mm
Scoop Height:	50mm
Nominal Weight:	580g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 28-Jul-20

Note: All dimensions are nominal