

**1. Part Number:** A336-1500

**2. Photograph:**



**3. Product Information:**

Description:	Open PharmaScoop® 1500ml
Material of Construction:	316L Stainless Steel
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	1500ml
Overall Length:	375mm
Body Length:	255mm
Body Diameter:	150mm
Scoop Height:	63mm
Nominal Weight:	580g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 28-Jul-20

Note: All dimensions are nominal