

1. Part Number: A336-2000

2. Photograph:



3. Product Information:

Description	Open PharmaScoop® 2000ml
Nominal Volume:	2000ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	420mm
Body Length:	295mm
Body Width:	165mm
Scoop Height:	63mm
Nominal Weight:	710g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 28-Jul-20

Note: All dimensions are nominal