


1. Part Number: 8640H-05S

2. Photograph:



3. Product Information:

| | |
|--|--|
| Description: | SteriWare® Volumetric Spoon – 5ml - Sterile |
| Use: | Single use, disposable measuring spoon |
| Material of Construction: | High Density Polyethylene (HDPE) – rigid, natural <ul style="list-style-type: none"> • FDA acceptable CFR 177.1520 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free |
| Method of Construction: | Single piece, injection moulded |
| Nominal Bowl Volume: | 5ml |
| Length: | 165mm |
| Bowl Diameter: | 29mm |
| Nominal Weight: | 4g |
| Moulding & Primary Packing Environment: | Class 100,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | Gamma Irradiation (25-45kGy) |
| Quantity per Box: | 100 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |
| ATEX: | This product has not been ATEX tested |

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 2-Mar-20