


1. Part Number: 5035T-1001

2. Photo:



3. Product Information:

Description	Liquid Master PTFE – 1000mm
Materials of Construction:	PTFE & FEP (for product contact parts) <ul style="list-style-type: none"> ○ FDA 21 CFR 177.1550 ○ USP Class VI ○ EU regulation 10/2011 ○ EC regulation 1935/2004
Outer Diameter of Outer Tube:	25mm
Internal Diameter of Outer Tube:	19mm
Max. Diameter of Tip:	32mm
Diameter of hole in end of tip:	10mm
Overall Length:	1000mm
BSE/TSE Free:	Yes
Nominal Weight:	585g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Oct-19