


1. Part Number: A416-10

2. Photograph:



Dimensions

A = 119mm

B = 98mm

C = 363mm

D = 212mm

Wall Thickness
= 1.2mm

3. Product Information:

| | |
|---------------------------------|--|
| Description: | Stainless Container 10L |
| Materials of Construction: | - Body: 316L stainless steel - Lid: 316L stainless steel - Clamp: 304 stainless steel - Gasket: Silicone (FDA acceptable grade) |
| Method of Construction: | All welds ground and polished Crevice free interior |
| Surface Finish: | Better than 0.5 microns Ra |
| Nominal Volume | 10 Litres |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this product are of vegetable origin. |
| Weight: | 3.6 Kg (body + lid + gasket + clamp) |
| Recommended Storage Conditions: | Dry and ambient temperature |

Note: The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Aug-20



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com