


1. Part Numbers: 5030A-1500

2. Photograph:



3. Product Information:

Description	Cup Sampler 100ml - 1500mm												
Volume of Cup	100ml (nominal)												
Material of Construction:	<table border="0"> <tr> <td>Cup</td> <td>316L stainless steel</td> </tr> <tr> <td>Cup Lid</td> <td>316L stainless steel</td> </tr> <tr> <td>Inner Rod</td> <td>316L stainless steel</td> </tr> <tr> <td>Outer Rod</td> <td>316L stainless steel</td> </tr> <tr> <td>Handles</td> <td>316L stainless steel</td> </tr> <tr> <td>Spring</td> <td>Stainless steel (not a product contact part)</td> </tr> </table>	Cup	316L stainless steel	Cup Lid	316L stainless steel	Inner Rod	316L stainless steel	Outer Rod	316L stainless steel	Handles	316L stainless steel	Spring	Stainless steel (not a product contact part)
Cup	316L stainless steel												
Cup Lid	316L stainless steel												
Inner Rod	316L stainless steel												
Outer Rod	316L stainless steel												
Handles	316L stainless steel												
Spring	Stainless steel (not a product contact part)												
Surface Finish:	Satin finish (better than 0.5µm Ra)												
Overall Length:	1550mm												
Diameter of Cup:	55mm												
Weight:	1.3kg												
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products												
Recommended Storage Conditions:	Dry and ambient temperature												

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20