


1. Part Number: 8055H-1000S

2. Photograph:



3. Product Information:

| | |
|---|--|
| Description: | SteriWare® ViscoThief - 1000mm |
| Use: | Single use, disposable liquid sampler (for high viscosity liquids) |
| Material of Construction: | High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • Conforms to REACH regulation 1907/2006 • BSE/TSE free O-ring Seal – virgin Silicone Rubber <ul style="list-style-type: none"> • Conforms to FDA CFR 177.2600 |
| Method of Construction: | Ultrasonic welding |
| Nominal Length | |
| Fully closed: | 991mm |
| Fully extended: | 1818mm |
| Diameter of Outer Tube: | 21mm |
| Diameter of hole in tip: | 10mm |
| Width of Handle: | 103mm |
| Nominal Weight of Sampler: | 158g |
| Assembly & Primary Packing Environment: | Class 100,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | Gamma Irradiation (25kGy minimum) |
| Number of Samplers per Box: | 20 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21