

**1. Part Number:** 8045P-125S

**2. Photograph:**



**3. Product Information:**


Description:	SteriWare PharmaScoop® PP 125ml - Sterile
Use:	Single use, disposable scoop
Material of Construction:	Polypropylene – rigid, opaque <ul style="list-style-type: none"> <li>• FDA acceptable CFR 177.1520</li> <li>• EU Regulations 10/2011 compliant</li> <li>• EC regulation 1935/2004 compliant</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Nominal Scoop Volume:	125ml
Overall Length:	196mm
Handle Length:	104mm
Maximum Width:	61mm
Nominal Weight:	27g
Moulding Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Scoops per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)



**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK.

Tel: +44(0)1675 466992 Email: [info@sampling.com](mailto:info@sampling.com) Web: [www.sampling.com](http://www.sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Dec-19