

1. Part Number: 8045H-1000S

2. Photograph:



3. Product Information:

Description:	SteriWare PharmaScoop® HDPE 1000ml - Sterile
Use:	Single use, disposable scoop
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • FDA acceptable CFR 177.1520 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Scoop Volume:	1000ml
Overall Length:	324mm
Handle Length:	149mm
Maximum Width:	150mm
Nominal Weight:	138g
Moulding Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Scoops per Box:	100 (supplied as 4 boxes of 25 scoops)
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21