


**1. Part Number:** 8015V-170S

**2. Photograph:**



### 3. Product Information:

|   |   |
|---|---|
| Description:                            | SteriWare® Very Narrow Micro Spatula 0.3ml - Sterile  |
| Use:                                    | Single use, powder spatula  |
| Material of Construction:               | High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1640</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul> |
| Method of Construction:                 | Single piece, injection moulded   |
| Nominal Volume:                         | 0.3ml   |
| Overall Length                          | 170mm   |
| Blade Length:                           | 40mm  |
| Blade Width:                            | 4mm   |
| Nominal Weight:                         | 3g  |
| Moulding & Primary Packing Environment: | Class 100,000 Cleanroom   |
| Individually Bagged?                    | Yes (heat sealed PE/PET bag)  |
| Method of Sterilisation:                | Gamma Irradiation (25 to 45 kGy)  |
| Number of Samplers per Box:             | 100   |
| Recommended Storage Conditions:         | Dry and ambient temperature   |
| Shelf Life:                             | 5 years from the date of manufacture  |
| ATEX:                                   | This device has not been ATEX tested  |

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Mar-20