


**1. Part Number:** 8015N-170S

**2. Photograph:**



**3. Product Information:**

Description:	SteriWare® Micro Spatula HDPE - Sterile
Use:	Single use, disposable spatula
Material of Construction:	High Impact Polystyrene (HIPS) <ul style="list-style-type: none"> <li>- FDA acceptable CFR 177.1640</li> <li>- EU Regulations 10/2011 compliant</li> <li>- EC Regulations 1935/2004 compliant</li> <li>- BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Nominal Volume:	0.5ml
Overall Length:	170mm
Spatula Blade Length:	40mm
Blade Width:	6mm
Nominal Weight:	4g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Scoops per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Mar-20