

# Sampler Specification

Specification Number

1030B

Revision

04



Part No.

1030B all lengths

Description

PowderThief - 19mm diameter

Part No:	Overall Length	Weight
1030B-600	680mm	1075g
1030B-1000	1080mm	1580g
1030B-1500	1580mm	2080g
1030B-2000	2080mm	2620g



Materials of Construction:

Inner Rod:	316L stainless steel
Handles:	316L stainless steel
Handle Boss:	316L stainless steel
Outer Tube:	316L stainless steel
Clutch Housing:	316L stainless steel
Clutch Socket:	316L stainless steel
Clutch:	PTFE
Disk:	PTFE

External Polish:

&lt;1 microns Ra

Internal Polish:

None

Recommended Storage Conditions:

Dry and ambient temperature

Spare Parts:

PTFE Clutch	Part No: SS002301
Clutch Housing	Part No: SS001901
PTFE Disk	Part No: SS002001

Specification Approved By:

Signature:

Name: Kuldeep Flora

Position: Quality Manager

Date: 17-Feb-20