


**1. Part Number:** 1030BT all sizes

**2. Product Information:**



Part No:	Nominal Volume	Weight
1030BT-2	2ml	173g
1030BT-3	3ml	163g
1030BT-5	5ml	146g
1030BT-8	8ml	124g
1030BT-10	10ml	112g

Description

Tips for use with 19mm diameter PowderThief Bodies

Materials of Construction:

316L stainless steel

Finish:

<1 microns Ra

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin

Recommended Storage Conditions:

Dry and ambient temperature

**3. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jan-21