

# Sampler Specification

Specification Number

1030C



Revision

04

Part No. 1030C all lengths

Description PowderThief - 25mm diameter

Part No:	Overall Length	Weight
1030C-600	680mm	1250g
1030C-1000	1080mm	1780g
1030C-1500	1580mm	2510g


**Materials of Construction:**

Inner Rod:	316 stainless steel
Handles:	316 stainless steel
Handle Boss:	316 stainless steel
Outer Tube:	316 stainless steel
Clutch Housing:	316 stainless steel
Clutch Socket:	316 stainless steel
Clutch:	PTFE
Disk:	PTFE

External Polish: &lt;1 microns Ra

Internal Polish: None

Recommended Storage Conditions: Dry and ambient temperature

Spare Parts:	PTFE Clutch	Part No: SS002301
	Clutch Housing	Part No: SS001902
	PTFE Disk	Part No: SS002002

Specification Approved By:	Signature: 	Name: Kuldeep Flora Position: Quality Manager Date: 17-Feb-20
Page 1 of 1	<b>Sampling Systems Ltd.</b> 4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com	