

Sampler Specification

Specification Number

1030D



Revision

05

Part No. 1030D all lengths

Description PowderThief - 32mm diameter

Part No:	Overall Length	Weight
1030D-600	680mm	1.6 Kg
1030D-1000	1080mm	2.3 Kg
1030D-1500	1580mm	3.2 Kg
1030D-2000	2080mm	4.0 Kg


Materials of Construction:

Inner Rod:	316L stainless steel
Handles:	316L stainless steel
Handle Boss:	316L stainless steel
Outer Tube:	316L stainless steel
Clutch Housing:	316L stainless steel
Clutch Socket:	316L stainless steel
Clutch:	PTFE
Disk:	PTFE

External Polish: <1 microns Ra

Internal Polish: None

Recommended Storage Conditions: Dry and ambient temperature

Spare Parts:	PTFE Clutch	Part No: SS002301
	Clutch Housing	Part No: SS001903
	PTFE Disk Ø32	Part No: SS002003

Specification Approved By:	Signature: 	Name: Kuldeep Flora Position: Quality Manager Date: 17-Feb-20
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com	