

1. Part Number: 8220P-03

2. Label Design:



3. Product Information:

Description: PharmaLabel - Quality Control – Yellow
 Adhesive: PharmaLabel Adhesive
 Colours: Yellow printing on a white label
 FDA Acceptable Adhesive? FDA CFR 175.105

The adhesive complies with the European food directives and legislations, FDA 175.105 and the German recommendations XXI as published by BfR. BfR (Bundesinstitut für Risikobewertung) is the German Federal Institute for Risk Assessment. The adhesive can be used in direct contact with dry and moist, non fatty foodstuffs.

The adhesive conforms to: European food regulations 1935/2004/EC, 10/2011/EU and DIN ISO 3826-1.

Dimensions: 95 x 95mm
 Number of Labels per Roll: 500 (approx).
 Temperature: Working Temperature Range -40 to +70°C (-40 to +160°F)
 Weight: 1010g per roll
 Method of Sterilisation: N/A
 Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 26-Sep-19



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
 Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com