

1. Part Number: 8220P-003

2. Label Design:



3. Product Information:

Description: PharmaLabel - Quality Control – Yellow

Adhesive: PharmaLabel Adhesive

Colours: Yellow printing on a white label

FDA Acceptable Adhesive? FDA CFR 175.105

The adhesive complies with the European food directives and legislations, FDA 175.105 and the German recommendations XXI as published by BfR. BfR (Bundesinstitut für Risikobewertung) is the German Federal Institute for Risk Assessment. The adhesive can be used in direct contact with dry and moist, non fatty foodstuffs.

The adhesive conforms to: European food regulations 1935/2004/EC, 10/2011/EU and DIN ISO 3826-1.

Dimensions: 150 x 150mm

Number of Labels per Roll: 250 (approx).



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

Temperature	Working Temperature Range -40 to +70°C (-40 to +160°F)
Weight:	1300g per roll
Method of Sterilisation:	N/A
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Sep-19