

1. Part Number: A426-50**2. Photograph:****3. Product Information:**

Description:	Mini-Pot 50ml
Method of Construction:	Crevice free
Material of Construction:	316L Stainless Steel
Nominal Volume:	50ml
Brimful Volume:	60ml
Overall Height (body & lid):	42mm
Body Height:	40mm
Maximum Body Diameter:	50mm
Lid Diameter:	53mm
Nominal Weight:	170g
	Bright Polish (Ra < 0.8 microns)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: If you require container with a watertight seal please see our containers type A836

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 19-Mar-21