

**1. Part Number:** 5005C-500

**2. Photo:**



### 3. Product Information:

|                             |   |
|-----------------------------|---|
| Part No.                    | 5005C-500   |
| Description                 | Stainless Steel Dipper – 500ml  |
| Material of Construction:   | 316L stainless steel  |
| Method of Construction:     | Fully Welded Construction   |
| Surface Finish:             | Cup better than 0.5 microns Ra<br>Handle better than 1.0 microns Ra                     |
| Nominal Cup Volume:         | 500ml   |
| Brimful Cup Volume:         | 520ml   |
| Overall Length:             | 1500mm  |
| Max. Width (cup and handle) | 100mm   |
| Cup Diameter:               | 90mm  |
| Cup Height:                 | 110mm   |
| Nominal Weight:             | 1250g   |
| BSE/TSE statement:          | All polishing compounds used in the manufacture of this product are of vegetable origin |

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Nov-19