

**1. Part Number:** A200-100

**2. Photograph:**



### 3. Product Information:

Description:	Pallet Knife with Plastic Handle
Material of Construction:	- Blade: Stainless Steel - Handle: Black PVC
Method of Construction:	Two piece construction
Surface Finish:	Better than 0.5 microns Ra
Blade Length:	100mm
Blade Width:	20mm
Overall Length:	190mm
Nominal Weight:	25g
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Apr-20