

1. Part Number: A200-300

2. Photograph:



3. Product Information:

Description: Pallet Knife with Plastic Handle
Material of Construction: - Blade: Stainless Steel
- Handle: Black PVC
Method of Construction: Two piece construction
Surface Finish: Better than 0.5 microns Ra
Blade Length: 300mm
Blade Width: 35mm
Overall Length: 415mm
Nominal Weight: 121g
Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Apr-20