

**1. Part Number:** A201-100

**2. Photograph:**



**3. Product Information:**

Description:	Pallet Knife with Wooden Handle
Material of Construction:	- Blade: Stainless Steel - Handle: Rosewood
Method of Construction:	Two piece construction
Surface Finish:	Better than 0.5 microns Ra
Blade Length:	105mm
Blade Width:	18mm
Overall Length:	190mm
Nominal Weight:	28g
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Apr-20