

1. Part Numbers: A264-6

2. Photograph:



3. Product Information:

Description	Stainless Bucket 6L
Nominal Volume:	6L
Brimful Volume:	7.3L
Material of Construction:	304 stainless steel
Method of Construction:	Crevice free interior
Surface Finish:	Brush Finish - better than 1.0 microns Ra
Max. Diameter:	257mm
Bucket Height:	240mm
Graduations:	Yes (every 1L)
Skirt Band:	No
Nominal Weight:	980g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jun-20