

1. Part Numbers: A264-12

2. Photograph:



3. Product Information:

Description	Stainless Bucket 12L
Nominal Volume:	12L
Brimful Volume:	12.5L
Material of Construction:	304 stainless steel
Method of Construction:	Crevice free interior
Surface Finish:	Brush Finish - better than 1.0 microns Ra
Max. Diameter:	311mm
Bucket Height:	290mm
Graduations:	Yes (every 1L)
Skirt Band:	Yes
Nominal Weight:	1.4 kg
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jun-20