


1. Part Number: A736-500

2. Photograph:



3. Product Information:

Description:	Stainless Beaker 500ml
Method of Construction:	Crevice free, all welds ground and polished
Material of Construction:	316L Stainless Steel
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	500ml
Dimensions:	
A - Internal Diameter:	82mm
B - Nominal Overall Diameter:	95mm
C - Nominal Height:	100mm
Wall Thickness	1.2mm
Nominal Weight:	275g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 30-Jun-21