

	<h1>Product Specification</h1>	Specification Number	8095B-150S
		Revision	01

Part No. 8095B-150S

Description Disposable Powder Spatula Blue – 150mm

Material of Construction: Polypropylene (PP) – rigid, blue

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulation 10/2011
- Conforms to EC Regulation 1935/2004



Method of Construction: Single piece, injection moulded

Nominal Length: 230mm

Blade Length: 145mm

Blade Width: 18mm

Nominal Weight: 16g

Moulding & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE bag)


Method of Sterilisation: Gamma Irradiation (25-45kGy)

Number of Samplers per Box: 100

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare is a registered trade marks belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 20-Jul-15
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	