

| | | | |
|-----------------------------------------------------------------------------------|--------------------------------|----------------------|------------|
|  | <h1>Product Specification</h1> | Specification Number | 8095B-250S |
| | | Revision | 01 |

Part No. 8095B-250S

Description Disposable Powder Spatula Blue – 250mm

Material of Construction: Polypropylene (PP) – rigid, blue

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulation 10/2011
- Conforms to EC Regulation 1935/2004



Method of Construction: Single piece, injection moulded

Nominal Length: 330mm

Blade Length: 246mm

Blade Width: 18mm

Nominal Weight: 22g

Moulding & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE bag)


Method of Sterilisation: Gamma Irradiation (25-45kGy)

Number of Samplers per Box: 100

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare is a registered trade marks belonging to Sampling Systems Ltd.

| | | |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|
| Specification Approved By: | Signature:  | Name: Jonathan Fenner Position: Director Date: 21-Apr-17 |
| Page 1 of 1 | Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill, Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com | |