

**1. Part Number:** 8890P-500S

**2. Photograph:**



**3. Product Information:**

Description:	SteriWare® Powder Bottle 500ml - Sterile
Use:	Single use, disposable bottle
Material of Construction:	Polypropylene (PP) – opaque <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulation 10/2011</li> <li>• Conforms to EC Regulation 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Nominal Volume:	500ml
Brimful Volume:	540ml
Method of Construction:	Single piece, injection moulded
Height:	106mm
Lid Diameter:	94mm
Nominal Weight:	83g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Bottles per Box:	37
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

ATEX:

This product has not been ATEX tested.

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-May-19