

1. Part Number: 8432P-500S

2. Photograph:



3. Product Information:

Description:	SteriWare® Graduated Jug 500ml - Sterile
Use:	Single use, disposable jug
Material of Construction:	Polypropylene (PP) – virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	500ml
Brimful Volume:	635ml
Top Diameter (excl. spout):	90mm
Max Width (inc. spout and handle)	135mm
Base Diameter:	75mm
Height:	130mm
Nominal Weight:	58g
Graduation Marks:	To 500ml (every 10ml)
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 - 45 kGy)
Quantity per Box:	35

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com**5. Country of Origin:**

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Mar-20