

	<h1>Product Specification</h1>	Specification Number	8209H-20S
		Revision	01

Part No. 8209H-20S

Description SteriWare Long Handled Spoon (Straight)

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004



Method of Construction: Single piece, injection moulded

Bowl Volume: 20ml

Handle Length: 285mm

Bowl Length: 88mm

Bowl Width: 38mm

Nominal Weight: 22g

Moulding & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)


Method of Sterilisation: Gamma Irradiation (min dose 25kGy)

Number of Samplers per Box: 100

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare is a registered trade marks belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 13-Jul-18
----------------------------	--	--