


**1. Part Number:** 8115N-150S

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® V-Spatula Type 2 - 150mm - Sterile
Use:	Single use, powder spatula
Material of Construction:	Polystyrene - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1640</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Blade Shape:	V
Blade Style:	Closed Tip
Overall Length	248mm
Blade Length:	147mm
Blade Width:	10mm
Nominal Weight:	11g
Nominal Volume (flat-filled):	2.3ml
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 45 kGy)
Number of Samplers per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

### 4. Change Control:

Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Oct-20