


1. Part Number: 8115A-250S

2. Photograph:



3. Product Information:

Description:	SteriWare® V-Spatula Type 2 - 250mm - Sterile
Use:	Single use, powder spatula
Material of Construction:	Polystyrene - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1640 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Blade Shape:	V
Blade Style:	Closed Tip
Overall Length	348mm
Blade Length:	247mm
Blade Width:	20mm
Nominal Weight:	19g
Nominal Volume (flat-filled):	25ml
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 45 kGy)
Number of Samplers per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Oct-20