

**1. Part Number:** M1000B-2500

**2. Photograph:**



**3. Product Information:**

Description:

Metal Detectable Scoop 2500ml - Blue

Material of Construction:

Polypropylene (PP)

- Metal detectable grade
- Only virgin material is used
- PP conforms to FDA 21 CFR 177.1520
- Blue colourant conforms to FDA 21 CFR 178.3297
- Filler conforms to FDA 21 CFR 186.1374
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Injection moulded as a single piece

Overall Length:

357mm

Handle Length:

147mm

Scoop Width:

166mm

Nominal Weight:

313g

Moulding Environment:

Class 100,000 Cleanroom

Packing Environment:

Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Number of scoops per Pack:

4

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare & Metal Detectable ranges please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 19-Apr-19