


1. Part Number: M1010B-250

2. Photograph:



3. Product Information:

Description:	Metal Detectable V-Spatula 250mm (Type1) - Blue
Material of Construction:	<p>Polypropylene (PP)</p> <ul style="list-style-type: none"> • Metal detectable grade • Only virgin material is used • PP conforms to FDA 21 CFR 177.1520 • Blue colourant conforms to FDA 21 CFR 178.3297 • Filler conforms to FDA 21 CFR 186.1374 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Injection moulded as a single piece
Overall Length:	345mm
Blade Length:	246mm
Blade Width:	20mm
Nominal Weight of Sampler:	28g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Number of samplers per Pack:	4
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare & Metal Detectable ranges please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Apr-19