

1. Part Number: M1023B-100

2. Photograph:



3. Product Information:

Description:

Metal Detectable Liquid Funnel 100mm - Blue

Material of Construction:

Polypropylene (PP)

- Metal detectable grade
- Only virgin material is used
- PP conforms to FDA 21 CFR 177.1520
- Blue colourant conforms to FDA 21 CFR 178.3297
- Filler conforms to FDA 21 CFR 186.1374
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Injection moulded as a single piece

Maximum Diameter:

117mm

Overall Height:

105mm

OD at Bottom of Spout

10mm

ID at Bottom of Spout:

8mm

OD at Top of Spout

17mm

Nominal Weight:

27g

Moulding Environment:

Class 100,000 Cleanroom

Packing Environment:

Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Number of funnels per Pack: 4
Recommended Storage Conditions: Dry and ambient temperature
Shelf Life: 5 years from the date of manufacture
ATEX: This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.
To be kept informed of all major changes to any product in the SteriWare & Metal Detectable ranges please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-May-19