

**1. Part Number:** 1610A-1500

**2. Photo:**



### 3. Product Information:

Description	Aluminium Slot Sampler 1500mm
Material of Construction:	Anodised Aluminium
Overall Length:	1500mm
Max. Sampling Depth:	1025mm
Number of Slots:	7
Slot Length:	100mm
Max. Slot Width:	13mm
Body Diameter:	25mm
Max. Sample Volume:	190ml
Nominal Weight:	1015g
Storage Conditions:	Room temperature

Note: The Tip is riveted to the outer tube and can not be removed

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 15-Feb-21