

# IKA

designed for scientists



## VX 2 Test tube attachment

/// Data Sheet

For up to 36 test tubes or centrifugal tubes with a diameter of 16 mm.

Dimensions (W x D x H): 140 x 145 x 115 mm

Material: macrolon

Weight: 300 g

[www.ika.com](http://www.ika.com)

Subject to technical changes





designed for scientists

### Technical Data

Width [mm]	87.5
Depth [mm]	30
Weight [kg]	0.309
Permissible ambient temperature [°C]	5 - 40