



designed for scientists



## AS 501.6 Separating funnel attachment

/// Data Sheet

For shaking out, eluting, extracting, gassing out, dissolving, enriching, etc. Adjustment for the clamping rolls is infinitely variable, the set-up height can be changed by means of clamping devices. This attachment will hold 4 x 1,000 ml separating funnels.

Dimensions (W x D x H): 480 x 505 x 225 mm

Weight: 5500 g



designed for scientists

### Technical Data

Width [mm]	84.6
Depth [mm]	56
Weight [kg]	5.588
Permissible ambient temperature [°C]	5 - 40