



## AS 4000.2 Fixing clip attachment

### /// Data Sheet

For shaking of flasks, Erlenmeyer flasks and bottles with a round crosssection (without fixing clips).

Dimensions (W x D x H): 470 x 440 x 25 mm

Number of fixing clips (volume):

Erlenmeyer flasks: 55 x AS 2.1 (25 ml), 52 x AS 2.2 (50 ml), 25 x AS 2.3 (100 ml), 16 x AS 2.4 (250 ml), 11 x AS 2.5 (500 ml), 7 x AS 2.6 (1000 ml), 4 x AS 2.7 (2000 ml)

Volumetric flasks: 55 x AS 2.1 (25 ml), 55 x AS 2.2 (50 ml), 32 x AS 2.3 (100 ml), 21 x AS 2.4 (250 ml), 21 x AS 2.5 (500 ml), 12 x AS 2.6 (1000 ml), 5 x AS 2.7 (2000 ml)

Set-up plate: 430 x 430 mm

Weight: 2.650 g

Clamping screws for fixing the attachment are included in the scope of delivery of the KS 4000 incubator shaker.



designed for scientists

For more information on our shaker accessories, please refer to the 'Selection Guide - Shakers and attachments' in the section 'Downloads'.



designed for scientists

### Technical Data

Depth [mm]	150
Height [mm]	60
Weight [kg]	1.315
Permissible ambient temperature [°C]	5 - 40