



Dimensions=mm

I. Product Description

Catalog number	#22500
Product description	Specimen container, 7mL, Biju, flat base, LxΦ=50 x 17mm, transparent, material of polystyrene, White crew cap, no label, bag pack

II. General Applications

Manufactured from virgin ultra clear Polystyrene.
Leak-proof at -95kPa (14psi);
Unique leak tight cap ensures exceptional leak proof seal.

III. Product Features

Label	No
Graduations	No
Cap type	Flat Screw cap
Base shape	Flat base
Working Volume	7mL
Autoclavable (max. °C)	No
Sterile	ETO sterile (SAL=10 ⁻⁶)

IV.Product Dimensions

(For more details, please check up above Drawing List)

Container Outer Diameter(O.D.) 17mm
Container Height 50mm

V.Material & colours

Product material Polystyrene
Colour of product Transparent
Cap material HDPE (High-density Polyethylene)
Colour of cap White

VI.Purity & certification

Satisfies the requirement IATA, -95kPa (14psi) / BS5213:1975(1989)
Product category IVD
Certificates Certified CE (98/79/EC) / FDA

VII.Packaging

Minimum order qty. 700
Type of inner packaging Bulk pack, 700 pieces/case
Piece(s) / inner pack No inner pack
Piece(s) / case 700
Case volume 0.02cmb (31*25*25,cm)
G.W. per case 3.8 Kgs

VIII.Related photos

/