


19 July 2019

TECHNICAL DATA SHEET SLE7538

Date of issue: 19 July 2019

Product Code: SLE7538

Description




60ml container, with Printed Label, in polystyrene, with a plastic screw cap with integral moulded seal. Product supplied in 300's, bagged in 100's. CE Marked. Class I IVD.

Declaration of Conformity

This product comes under the remit of the In Vitro Medical Devices Directive 98/79/EC and is classed as a general IVD product under annex III (specimen receptacle as per article 1(b) of the directive).


As the manufacturer we hereby declare that the requirements and obligations, as described in annex III sections 2 to 5 of the Directive have been met.



QADOCPROD: SLE7538 PAGE 1 of 2 Authorised by: Andrew Littlewood

Microbiological Status	Manufactured using an automated aseptic manufacturing system under laminar flow HEPA filtration.
Standard to which the product conforms	BS EN 14254:2004
Material of construction (body)	Polystyrene
Material of construction (Cap shell)	High Density Polyethylene
Volume (ml)	60ml
Bottle Dimensions (mm)	64.6 x39.0mm (approx)
Operational temperature Range	-10 to 60 Degrees Centigrade
Inner Pack Quantity	100 in re sealable gripper bag
Outer Box Quantity	300

QADOCPROD: SLE7538 PAGE 2 of 2, ISSUE DATE 19.07.2019



Authorised by: Andrew Littlewood

