


19 July 2019

TECHNICAL DATA SHEET SLE7540

Date of issue: 19 July 2019

Product Code: SLE7540


Description



60ml container, in polystyrene, with a plastic screw cap with integral moulded seal. Plain label. Product supplied in 300's, bagged in 100's.

Declaration of Conformity

I hereby declare that the product described in this document has been manufactured to the requirements of BS EN14254:2004



Authorised by: Andrew Littlewood

QADOCPROD: SLE7540 PAGE 1 of 2

Microbiological Status	Manufactured using an aseptic manufacturing system under laminar flow HEPA filtration.
Standard to which the product conforms	BS EN 14254:2004
Material of construction (body)	Polystyrene
Material of construction (Cap)	High Density Polyethylene
Volume (ml)	60ml
Bottle Dimensions (mm)	64.6 x 39.0 mm (approx.)
Operational temperature Range	-10 to 60 Degrees Centigrade
Inner Pack Quantity	100 in re sealable gripper bag
Outer Box Quantity	300

QADOCPROD: SLE7540 PAGE 2 of 2, ISSUE DATE: 19 July 2019

Authorised by: Andrew Littlewood

