

18 July 2019

## TECHNICAL DATA SHEET SLE7554

Date of issue: 18 July 2019

Product Code: SLE7554

### Description




100ml container, in polystyrene, with a metal screw cap with integral flowed in seal. Plain label. Product supplied in 160's, bagged in 40's.

### Declaration of Conformity


I hereby declare that the product described in this document has been manufactured to the requirements of BS EN14254:2004

Authorised by: Andrew Littlewood

A handwritten signature in black ink, appearing to be 'AL', written over a white background.

Microbiological Status	Manufactured using an aseptic manufacturing system under laminar flow HEPA filtration.
Standard to which the product conforms	BS EN 14254:2004
Material of construction (body)	Polystyrene
Cap Liner Material	Daraseal
Material of construction (Cap shell)	Tinplate
Volume (ml)	100ml
Bottle Dimensions (mm)	77 x 49mm (approx.)
Operational temperature Range	-10 to 60 Degrees Centigrade
Inner Pack Quantity	40 in re sealable gripper bag
Outer Box Quantity	160

QADOCPROD: SLE7554 PAGE 2 of 2, ISSUE DATE: 18 July 2019



Authorised by: Andrew Littlewood

