

Material description: ACETIC ACID GLACIAL A.R.

CAS No: 64-19-7

<i>Supplied as:</i>	Product code	Description
	CHE1016	ACETIC ACID GLACIAL A.R. 2.5L.
	CHE1018	ACETIC ACID GLACIAL A.R. 500ml.
	CHE1020	ACETIC ACID GLACIAL A.R. 25L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.8			
Boiling range	°C	117.5	119.0		
Wt. per ml.	g	1.048	1.051		
Chloride	%		0.0001		
Sulphate	%		0.0001		
Iron	%		0.00005		
Heavy Metals; as Pb	%		0.00005		
Acetaldehyde	%		0.005		
Formic Acid	%		0.05		
Water	% w/w		0.2		

Revision date: 10/08/2004

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0011/MS/20040810