

Material description: ACETONE A.R.

CAS No: 67-64-1

<i>Supplied as:</i>	Product code	Description
	CHE1034	ACETONE A.R. 2.5L.
	CHE1036	ACETONE A.R. 25L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Appearance					Clear colourless, free from suspended material
Purity	% w/w	99.8			
Wt. per ml.	g	0.789	0.793		
Colour	Hazen		5		
Water	% w/w		0.3		
Distillation IBP	°C	55.6			
Distillation DP	°C		56.6		

Revision date: 28/07/2023

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0015/MS/20230728