

Material description: ACETONITRILE S HPLC (far uv grade)

CAS No: 75-05-8

*Supplied as:* Product code Description  
CHE1054 ACETONITRILE S HPLC (far uv grade) 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.9			
Acidity; acetic	%		0.005		
Free alkali	%		0.001		
Water	% w/w		0.02		
Non-volatiles	% w/v		0.0002		
Transmission; 195nm	%	80			
Transmission; 210nm	%	90			
Transmission; 240nm	%	99			
Transmission; 254nm	%	95			

Revision date: 10/06/1999

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2227/MS/19990610