



Material description: ALKALINE PYROGALLOL (Soln B) (Potassium Hydroxide)
CAS No: 1310-58-3

Supplied as: Product code Description
CHE1088 ALKALINE PYROGALLOL (Soln B) (Potassium Hydroxide) 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Potassium hydroxide content	mol/L	3.544	3.616	3.580	

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3114/MS/19890101