



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: 1-AMINOBTANE pure

CAS No: 109-73-9

Supplied as: Product code Description
CHE1120 1-AMINOBTANE pure 250ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Boiling range	°C	76	78		
Wt. per ml.	g	0.736	0.738		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1667/MS/19890101