



Material description: AMMONIA SOLUTION 25% w/w
CAS No: 1336-21-6

Supplied as: Product code Description
CHE1124 AMMONIA SOLUTION 25% w/w 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	% w/w	23	26		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2970/MS/19890101